

AIMScribe v3 — Where Everything Lives

Access sheet. URLs, field names, and the file that holds each value. No passwords are printed here.

Every value named below is in **deploy/local/.env** on the bench machine. The UIU server has its own **deploy/uiu/.env** with different values, generated there. Neither file is in git — both AIMScribe repositories are public.

1 · Operations dashboard

What	Value / where
URL	http://localhost:6060/api/v2/dashboard
Asks for	the administrator key
Key lives in	.env → AIMS_ADMIN_KEY

The key is held in memory only, never in browser storage, so it is re-entered after a reload.

2 · Database viewer (Adminer)

Field	What to enter
URL	http://localhost:8080
System	PostgreSQL
Server	postgres
Username	aims_admin
Password	.env → POSTGRES_SUPERUSER_PASSWORD
Database	aims_recordings or aims_clinical

The Server field must be postgres, never localhost. Adminer runs inside a container, so localhost there is Adminer itself: the connection is refused and the page returns empty with no useful error. Start it with `docker compose --profile db up -d dbviewer`. It listens on the server machine only.

3 · PostgreSQL, direct

Role	Username / password variable
Host, port	localhost : 5433
Superuser	aims_admin / POSTGRES_SUPERUSER_PASSWORD
Recordings	aims_recordings / AIMS_RECORDINGS_PASSWORD
Clinical	aims_clinical_writer / AIMS_CLINICAL_PASSWORD
Monitor, read-only	aims_monitor / AIMS_MONITOR_PASSWORD

Two databases, two roles, and neither can reach the other's (SRS-DBA-20). **aims_recordings** holds sessions, integrity chains and receipts, and contains no patient names by design. **aims_clinical** holds the patient data.

4 · Command line — nothing to type wrong

Command	Shows
python tools/db.py tables	every table in both databases, real counts
python tools/db.py sessions	the most recent recordings
python tools/db.py patient <id>	one patient across both databases

python tools/db.py sql "select ..."	any read; --db clinical for the other
python tools/db.py analyze	refresh row estimates after a crash

It goes through the database container directly, so there is no hostname to get wrong and no password to enter. Reads only, unless --write is passed.

5 · Storage, queue and the recorder

Service	Address / where the secret is
Object storage console	http://localhost:9001
Storage user, secret	MINIO_ACCESS_KEY, MINIO_SECRET_KEY
Redis	localhost : 6379 / REDIS_PASSWORD
Recorder, doctor's PC	ws://127.0.0.1:5050/ws (loopback only)
Recorder local API key	that PC's .env → AIMS_LOCAL_API_KEY

6 · Integration keys

Key	Where it comes from
CMED Channel B	.env → AIMS_CMED_KEY · mint: POST /api/v2/admin/cmed-key
Archive worker	.env → AIMS_WORKER_KEY
Enrolment token	deploy/local/out/recorder/enrollment.token
Mint a token	POST /api/v2/admin/enrollment-token (X-Admin-Key)

One enrolment token binds one laptop to one clinic, once. Fourteen rooms means fourteen tokens. A CMED key is shown only at the moment it is minted.

7 · Ports

Port	Service
6060	AIMS LAB server — API and dashboard
8080	Adminer, server machine only, --profile db
5433	PostgreSQL
9001	Object storage console
5050	The recorder on a doctor's PC, loopback only
3000	CMED test site, only with --cmed

Not port 6000. Browsers and Node both refuse it — it belongs to X11. The server uses 6060 inside and 443 at the edge.

8 · The archive on disk

What	Path
Archive root	AIMS_ARCHIVE_PATH (/srv/aims/archive on the UIU server)
Databases	AIMS_DB_PATH (/srv/aims/db)
Backups	AIMS_BACKUP_PATH (/srv/aims/backups)
One consultation	<clinic>/<doctor>/<date>/<name>/

Inside that folder	the .wav, its .manifest.json, its clinical .json
---------------------------	--

The dashboard's recordings table shows this folder per recording with a copy button, so removing one by hand is a copy and a paste rather than a query written out.