

AIMScribe v3 — Where Everything Lives

Access sheet. URLs, field names, and the file that holds each value. No passwords are printed here.

Every value named below is in **deploy/local/.env** on the bench machine. The UIU server has its own **deploy/uiu/.env** with different values, generated there. Neither file is in git — both AIMScribe repositories are public.

1 · Every page, in one place

Page	Address	Needs
Operations dashboard	http://localhost:6060/api/v2/dashboard	the administrator key
CMED dummy page	http://localhost:3000	nothing — this is the page a doctor uses
CMED protocol bench	http://localhost:3000/protocol-test	nothing; shows every message and reply
Database viewer	http://localhost:8080	the superuser, and Server = postgres
Object storage console	http://localhost:9001	the storage user
Server health	http://localhost:6060/health	nothing
The recorder, on a doctor's PC	ws://127.0.0.1:5050/ws	loopback only; the page connects, not a person

The CMED page only runs when the stack was started with `--cmcd`. It is the one to open to record a real consultation: open a patient, speak, press Build Prescription, then open the next patient or Stop.

The dashboard's administrator key is held in memory only, never in browser storage, so it is re-entered after a reload.

2 · Database viewer (Adminer)

Field	What to enter
URL	http://localhost:8080
System	PostgreSQL
Server	postgres
Username	aims_admin
Password	.env → POSTGRES_SUPERUSER_PASSWORD
Database	aims_recordings or aims_clinical

The Server field must be postgres, never localhost. Adminer runs inside a container, so localhost there is Adminer itself: the connection is refused and the page returns empty with no useful error. Start it with `docker compose --profile db up -d dbviewer`. It listens on the server machine only.

3 · The logins, one row each

Host and port for all three database roles: localhost : 5433. The value of each secret is in **deploy/local/.env** under the variable named here — copy it into the last column once, and this page is the whole set.

What it opens	Username	Variable holding the secret	Value
Dashboard	—	AIMS_ADMIN_KEY	

Adminer / superuser	aims_admin	POSTGRES_SUPERUSER_PASSWORD	
aims_recordings	aims_recordings	AIMS_RECORDINGS_PASSWORD	
aims_clinical	aims_clinical_writer	AIMS_CLINICAL_PASSWORD	
Monitor, read-only	aims_monitor	AIMS_MONITOR_PASSWORD	
Object storage	aimslocal	MINIO_SECRET_KEY	
Redis	—	REDIS_PASSWORD	
CMED Channel B	—	AIMS_CMED_KEY	
Archive worker	—	AIMS_WORKER_KEY	

To read them all at once on the bench machine:

```
cd deploy\local && findstr /R "PASSWORD KEY" .env
```

Two databases, two roles, and neither can reach the other's (SRS-DBA-20). **aims_recordings** holds sessions, integrity chains and receipts, and contains no patient names by design. **aims_clinical** holds the patient data.

4 · Command line — nothing to type wrong

Command	Shows
python tools/db.py tables	every table in both databases, real counts
python tools/db.py sessions	the most recent recordings
python tools/db.py patient <id>	one patient across both databases
python tools/db.py sql "select ..."	any read; --db clinical for the other
python tools/db.py analyze	refresh row estimates after a crash

It goes through the database container directly, so there is no hostname to get wrong and no password to enter. Reads only, unless --write is passed.

5 · Storage, queue and the recorder

Service	Address / where the secret is
Object storage console	http://localhost:9001
Storage user, secret	MINIO_ACCESS_KEY, MINIO_SECRET_KEY
Redis	localhost : 6379 / REDIS_PASSWORD
Recorder, doctor's PC	ws://127.0.0.1:5050/ws (loopback only)
Recorder local API key	that PC's .env → AIMS_LOCAL_API_KEY

6 · Integration keys

Key	Where it comes from
CMED Channel B	.env → AIMS_CMED_KEY · mint: POST /api/v2/admin/cmed-key
Archive worker	.env → AIMS_WORKER_KEY
Enrolment token	deploy/local/out/recorder/enrollment.token

Mint a token	POST /api/v2/admin/enrollment-token (X-Admin-Key)
---------------------	---

One enrolment token binds one laptop to one clinic, once. Fourteen rooms means fourteen tokens. A CMED key is shown only at the moment it is minted.

7 · Ports

Port	Service
6060	AIMS LAB server — API and dashboard
8080	Adminer, server machine only, --profile db
5433	PostgreSQL
9001	Object storage console
5050	The recorder on a doctor's PC, loopback only
3000	CMED test site, only with --cmed

Not port 6000. Browsers and Node both refuse it — it belongs to X11. The server uses 6060 inside and 443 at the edge.

8 · The archive on disk

What	Path
Archive root	AIMS_ARCHIVE_PATH (/srv/aims/archive on the UIU server)
Databases	AIMS_DB_PATH (/srv/aims/db)
Backups	AIMS_BACKUP_PATH (/srv/aims/backups)
One consultation	<clinic>/<doctor>/<date>/<name>/
Inside that folder	the .wav, its .manifest.json, its clinical .json

The dashboard's recordings table shows this folder per recording with a copy button, so removing one by hand is a copy and a paste rather than a query written out.